

Force plates

Exploration of semi-active control laws for the Two Axis aDaptable Ankle (TADA)

- The Two Axis aDaptable Ankle (TADA) is a semi-active prosthetic ankle that allows sagittal and frontal ankle angle control
- Method: walk at self-selected speed using TADA angles of neutral angle (90° dorsiflexion and no eversion), and from neutral. 5° dorsiflexion, 5° plantarflexion, 5° inversion, and 5° eversion (5 total angles)
- The TADA will have an IMU on the foot, hall sensors and encoders for the motors, and Europa+ for the load cell


Aim to investigate

- How to move the foot during the swing phase to increase minimum toe clearance?
- What is the right way for the TADA to match the slope?
- How to use ankle moment to select ankle angles?




Kieran Nichols Advisor: Peter Adamczyk Co-author: Becca Roembke





Toe-lift control



Anticipatory Slope-Matching Control

Moment-Limiting Control

